

For Research Use Only

EPHX2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86790-2-PBS



Basic Information

Catalog Number: 86790-2-PBS	GenBank Accession Number: BC013874	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2053	CloneNo.: 251753F5
Source: Rabbit	UNIPROT ID: P34913	
Isotype: IgG	Full Name: epoxide hydrolase 2, cytoplasmic	
Immunogen Catalog Number: AG1283	Calculated MW: 63 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86790-2-PBS targets EPHX2 as part of a matched antibody pair:

MP02802-1: 86790-3-PBS capture and 86790-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EPHX2 concentration was determined to be 15.72 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 24.32 ng/mL in HEK-293T cell extract based on a 1.3 mg/mL extract load and 11.98 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02802-1, Human EPHX2 Recombinant Matched Antibody Pair - PBS only. 86790-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag1283. 86790-2-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL