

For Research Use Only

LEPREL2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86813-2-PBS



Basic Information

Catalog Number: 86813-2-PBS	GenBank Accession Number: BC017217	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10536	CloneNo.: 251845C8
Source: Rabbit	UNIPROT ID: Q8IVL6	
Isotype: IgG	Full Name: leprecan-like 2	
Immunogen Catalog Number: AG8850	Calculated MW: 736 aa, 82 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86813-2-PBS targets LEPREL2 as part of a matched antibody pair:

MP03545-1: 86813-3-PBS capture and 86813-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean P3H3 concentration was determined to be 1,831.68 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load, 715.68 ng/mL in Ramos cell extract based on a 1.20 mg/mL extract load and 798.43 ng/mL in Raji cell extract based on a 1.50 mg/mL extract load.



Sandwich ELISA standard curve of MP03545-1, Human P3H3 Recombinant Matched Antibody Pair - PBS only. 86813-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8850. 86813-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.