

For Research Use Only

SH3PXD2A Recombinant monoclonal antibody

Catalog Number: 86832-2-RR



Basic Information

Catalog Number: 86832-2-RR	GenBank Accession Number: NM_001394015.1	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 9644	CloneNo.: 251694C3
Source: Rabbit	UNIPROT ID: Q5TCZ1	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: SH3 and PX domains 2A	
Immunogen Catalog Number: AG39231	Calculated MW: 125kDa 1133aa	
	Observed MW: 130-150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, NIH/3T3 cells, HSC-T6 cells
Species Specificity: human, mouse, rat	

Background Information

SH3PXD2A is also named as FISH, KIAA0418, SH3MD1, TKS5. SH3PXD2A contains an amino-terminal PX domain followed by five SH3 domains. It is a cytoplasmic protein in normal fibroblasts (PMID: 19464300). The p140 and p130 forms of SH3PXD2A may be generated by other splice variations. The complexity of SH3PXD2A isoforms is also evident in different cell types. For example in human platelets a single band of 150 kDa was detected, whereas all three forms were found in human fibroblasts and vascular smooth muscle cells (PMID: 9687503).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86832-2-RR (SH3PXD2A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86832-2-RR against Human SH3PXD2A were performed. The affinity constant is 0.906 nM.