

For Research Use Only

RAB34 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86839-3-PBS



Basic Information

Catalog Number: 86839-3-PBS	GenBank Accession Number: BC016841	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 83871	CloneNo.: 251867B8
Source: Rabbit	UNIPROT ID: Q9BZG1	
Isotype: IgG	Full Name: RAB34, member RAS oncogene family	
Immunogen Catalog Number: AG26696	Calculated MW: 251 aa, 28 kDa	
	Observed MW: 29 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

RAB34, a Rab protein family member, is involved in repositioning of lysosomes, activation of micropinocytosis, and protein transport. RAB34 is aberrantly expressed in various cancers and exhibits oncogenic properties. Recently, Rab34 is identified as a ciliary protein. Rab34 is required for internal assembly of cilia but not for cilia built on the cell surface.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86839-3-RR (RAB34 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86839-3-PBS in a different storage buffer formulation.