

For Research Use Only

SHANK1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86842-1-PBS



Basic Information

Catalog Number:

86842-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_016148

GeneID (NCBI):

50944

UNIPROT ID:

Q9Y566

Full Name:

SH3 and multiple ankyrin repeat domains 1

Calculated MW:

225 kDa

Observed MW:

159-225 kDa

Purification Method:

Protein A purification

CloneNo.:

251628B1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86842-1-RR (SHANK1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86842-1-PBS in a different storage buffer formulation.