

For Research Use Only

MMP10 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86845-1-PBS



Basic Information

Catalog Number: 86845-1-PBS	GenBank Accession Number: NM_002425.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4319	CloneNo.: 251067E10
Source: Rabbit	UNIPROT ID: P09238	
Isotype: IgG	Full Name: matrix metalloproteinase 10 (stromelysin 2)	
Immunogen Catalog Number: EG3038	Calculated MW: 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

86845-1-PBS targets MMP10 as part of a matched antibody pair:

MP02702-1: 86845-2-PBS capture and 86845-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HT-1080 human fibrosarcoma cells were cultured in DMEM with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were stimulated with PMA. An aliquot of the cell culture supernate was removed, assayed for human Total MMP-10, and measured 20.5 ng/mL



Cytometric bead array standard curve of MP02702-1, Total MMP-10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86845-2-PBS. Detection antibody: 86845-1-PBS. Standard: Eg3038. Range: 0.625-80 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 86845-1-PBS against Human MMP10 were performed. The affinity constant is 0.924 nM.