

For Research Use Only

HBZ Recombinant monoclonal antibody

Catalog Number: 86879-3-RR



Basic Information

Catalog Number: 86879-3-RR	GenBank Accession Number: BC027892	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3050	CloneNo.: 251873A7
Source: Rabbit	UNIPROT ID: P02008	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:200-1:800
Isotype: IgG	Full Name: hemoglobin, zeta	
Immunogen Catalog Number: AG11061	Calculated MW: 142 aa, 16 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : K-562 cells, IF/ICC : K-562 cells,
Species Specificity: human	

Background Information

HBZ, also named as Zeta-globin and HBAZ, belongs to the globin family. HBZ is an alpha-type chain of mammalian embryonic hemoglobin, synthesized primarily in the yolk sac. HBZ was first identified as an antisense viral protein that associated with cAMP-response element binding protein-2 (CREB-2) in HTLV-1-infected cells. HBZ is the only viral gene that is expressed in all ATL cases (PMID: 22787458). ShRNA repression of HBZ expression in established HTLV-1-transformed cell lines and newly immortalized T-lymphocytes also significantly suppressed T-lymphocyte proliferation (PMID: 19650892, 16917009).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86879-3-RR (HBZ antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed K-562 cells using HBZ antibody (86879-3-RR, Clone: 251873A7) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



Bi-layer interferometry (BLI) kinetic assays of 86879-3-RR against Human HBZ were performed. The affinity constant is below 1 pM.