

For Research Use Only

CCT2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86898-2-PBS



Basic Information

Catalog Number: 86898-2-PBS	GenBank Accession Number: BC113516	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10576	CloneNo.: 251618F8
Source: Rabbit	UNIPROT ID: P78371	
Isotype: IgG	Full Name: chaperonin containing TCP1, subunit 2 (beta)	
Immunogen Catalog Number: AG19264	Calculated MW: 535 aa, 57 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

86898-2-PBS targets CCT2 as part of a matched antibody pair:

MP02720-1: 86898-2-PBS capture and 86898-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CCT2 concentration was determined to be 600.11 ng/mL in K-562 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02720-1, Human CCT2 Recombinant Matched Antibody Pair - PBS only. 86898-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19264. 86898-1-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL