

For Research Use Only

HNRNPUL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86926-1-PBS



Basic Information

Catalog Number: 86926-1-PBS	GenBank Accession Number: BC002564	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 11100	CloneNo.: 252010C6
Source: Rabbit	UNIPROT ID: Q9BUJ2	
Isotype: IgG	Full Name: heterogeneous nuclear ribonucleoprotein U-like 1	
Immunogen Catalog Number: AG0899	Calculated MW: 120 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

86926-1-PBS targets HNRNPUL1 as part of a matched antibody pair:

MP02774-1: 86926-2-PBS capture and 86926-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HNRNPUL1 concentration was determined to be 87.9 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 95.9 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 33.9 ng/mL in HuH-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02774-1, HNRNPUL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86926-2-PBS. Detection antibody: 86926-1-PBS. Standard: Ag0899. Range: 0.312-40 ng/mL.