

For Research Use Only

HNRNPUL1 Recombinant monoclonal antibody

Catalog Number:86926-3-RR



Basic Information

Catalog Number: 86926-3-RR	GenBank Accession Number: BC002564	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 11100	CloneNo.: 252010C1
Source: Rabbit	UNIPROT ID: Q9BUJ2	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: heterogeneous nuclear ribonucleoprotein U-like 1	
Immunogen Catalog Number: AG0899	Calculated MW: 120 kDa	
	Observed MW: 120 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HeLa cells, MCF-7 cells, HEK-293 cells, mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Background Information

HNRNPUL1, also named as E1BAP5, HNRPUL1 and E1B-AP5, acts as a basic transcriptional regulator. It represses basic transcription driven by several virus and cellular promoters. HNRNPUL1 plays also a role in mRNA processing and transport. It is a cellular protein that interacts with the adenovirus protein E1B-55 and implicated in mRNA export. HNRNPUL1 is an hnRNP-like protein that likely mediates the interaction of NXF1 with mRNAs.(PMID:17267598). It may modulate PAN expression in complex with KSHV ORF57(PMID:22043172). HNRNPUL1 specifically bound to the Rae1-Nup98 complex independently of M protein, and this interaction was found to be RNA dependent, as it was greatly diminished in the presence of RNase.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86926-3-RR (HNRNPUL1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86926-3-RR against Human HNRNPUL1 were performed. The affinity constant is 0.707 nM.