

For Research Use Only

IL-16 Recombinant monoclonal antibody

Catalog Number: 86927-2-RR



Basic Information

Catalog Number: 86927-2-RR	GenBank Accession Number: AAC12732.1	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3603	CloneNo.: 251661G11
Source: Rabbit	Full Name: interleukin 16 (lymphocyte chemoattractant factor)	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Calculated MW: 9 kDa	
Immunogen Catalog Number: EG4831	Observed MW: 80 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Recombinant protein, Jurkat cells, HuT 78 cells
Species Specificity: human	

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86927-2-RR (IL-16 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86927-2-RR (IL-16 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.