

For Research Use Only

TRADD Recombinant monoclonal antibody

Catalog Number: 86932-1-RR



Basic Information

Catalog Number: 86932-1-RR	GenBank Accession Number: BC004491	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 8717	CloneNo.: 251976C2
Source: Rabbit	UNIPROT ID: Q15628	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: TNFRSF1A-associated via death domain	
Immunogen Catalog Number: AG7768	Calculated MW: 34 kDa	
	Observed MW: 34 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HeLa cells, SH-SY5Y cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 86932-1-RR (TRADD antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.