

For Research Use Only

# RCOR1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86954-3-PBS



## Basic Information

<b>Catalog Number:</b> 86954-3-PBS	<b>GenBank Accession Number:</b> NM_015156	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 23186	<b>CloneNo.:</b> 252115G5
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UKLO	
<b>Isotype:</b> IgG	<b>Full Name:</b> REST corepressor 1	
<b>Immunogen Catalog Number:</b> AG26770	<b>Calculated MW:</b> 53 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
human

## Product Information

86954-3-PBS targets RCOR1 as part of a matched antibody pair.

MP02924-1: 86954-3-PBS capture and 86954-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean RCOR1 concentration was determined to be 903.00 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load, 279.00 ng/mL in K-562 cell extract based on a 1.20 mg/mL extract load, 566.90 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 267.70 ng/mL in HuH-7 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02924-1, Human RCOR1 Recombinant Matched Antibody Pair - PBS only. 86954-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag26770. 86954-2-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.