

For Research Use Only

# CNTN3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86965-1-PBS



## Basic Information

Catalog Number:

86965-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5928

GenBank Accession Number:

NM\_020872.2

GeneID (NCBI):

5067

UNIPROT ID:

Q9P232

Full Name:

contactin 3 (plasmacytoma associated)

Calculated MW:

113 kDa

Observed MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

251948A2

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86965-1-RR (CNTN3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86965-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 86965-1-RR (CNTN3 antibody) at dilution of 1:500 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86965-1-PBS in a different storage buffer formulation.