

For Research Use Only

HP1BP3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86974-2-PBS



Basic Information

Catalog Number: 86974-2-PBS	GenBank Accession Number: BC045660	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 50809	CloneNo.: 252095F2
Source: Rabbit	UNIPROT ID: Q5SSJ5	
Isotype: IgG	Full Name: heterochromatin protein 1, binding protein 3	
Immunogen Catalog Number: AG18709	Calculated MW: 553 aa, 61 kDa	
	Observed MW: 75 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

HP1BP3 (heterochromatin protein 1-binding protein 3) is an evolutionarily conserved nuclear protein that functions as a non-canonical linker histone. It is broadly expressed in most vertebrate tissues and is structurally related to the histone H1 family, containing three globular domains and a highly basic C-terminal tail reminiscent of H1 proteins. Unique to HP1BP3 are an acidic N-terminal stretch and a consensus HP1-binding motif (PxVxL), through which it associates with heterochromatin protein 1 (HP1 α/β) and is enriched at perinucleolar heterochromatin regions marked by H3K9me3. Functionally, HP1BP3 acts as a modulator of gene expression and genome stability. Alternative splicing was reported to produce 4 different isoforms of HP1BP3, the calculated MW corresponding to 61, 57, 44, 15 kDa.

Storage

Storage:
Store at -80°C.


Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Various lysates were subjected to SDS PAGE followed by western blot with 86974-2-RR (HP1BP3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86974-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86974-2-RR against Human HP1BP3 were performed. The affinity constant is 4.61 nM.