

For Research Use Only

# ZBTB48 Recombinant monoclonal antibody

Catalog Number: 86981-6-RR



## Basic Information

<b>Catalog Number:</b> 86981-6-RR	<b>GenBank Accession Number:</b> BC013573	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul	<b>GeneID (NCBI):</b> 3104	<b>CloneNo.:</b> 252100E5
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P10074	<b>Recommended Dilutions:</b> WB: 1:1000-1:4000
<b>Isotype:</b> IgG	<b>Full Name:</b> zinc finger and BTB domain containing 48	
<b>Immunogen Catalog Number:</b> AG20301	<b>Calculated MW:</b> 688 aa, 77 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB: HEK-293 cells, HepG2 cells, Jurkat cells, K-562 cells
<b>Species Specificity:</b> human	

## Background Information

ZBTB48, also named as TZAP and HKR3, is a 688 amino acid protein, which contains 1 BTB (POZ) domain and 11 C2H2-type zinc fingers and belongs to the krueppel C2H2-type zinc-finger protein family. ZBTB48 is detected in adrenal gland and neuroblastoma. ZBTB48 binds to and regulates the J and/or S elements in MHC II promoter. Zinc finger and BTB domain containing 48 (ZBTB48), also known as telomeric zinc-finger associated protein (TZAP), is a protein that triggers the trimming of telomeres in the absence of shelterin. ZBTB48, is a relatively uncharacterized POK family protein with a POZ-domain and 11 zinc finger domains. ZBTB48, is ubiquitously expressed in human tissues. The ZBTB48 gene is mapped to human chromosome 1p36, which is commonly rearranged (leiomyoma and leukemias) or deleted in various cancers (neuroblastoma, melanoma, Merkel cell carcinomas, pheochromocytoma and breast and colon carcinomas). A correlation between deletions or rearrangements of HKR3 and human cancer suggest a role for ZBTB48 as a potential tumor suppressor.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
**Aliquoting is unnecessary for -20°C storage**

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86981-6-RR (ZBTB48 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86981-6-RR against Human ZBTB48 were performed. The affinity constant is 0.588 nM.