

TERF2 Recombinant monoclonal antibody

Catalog Number: 87020-1-RR

Basic Information

Catalog Number:

87020-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16798

GenBank Accession Number:

BC024890

GeneID (NCBI):

7014

UNIPROT ID:

Q15554

Full Name:

telomeric repeat binding factor 2

Calculated MW:

56 kDa

Observed MW:

65-69 kDa

Purification Method:

Protein A purification

CloneNo.:

252147A6

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, HeLa cells, K-562 cells, Daudi cells, MCF-7 cells, NIH/3T3 cells, HSC-T6 cells

Background Information

TRF2, also named as TRBF2 and Telomeric repeat-binding factor 2, is a 542 amino acid protein, which contains 1 HTH myb-type DNA-binding domain and localizes in the Nucleus. TRF2 binds the telomeric double-stranded 5'-TTAGGG-3' repeat and plays a central role in telomere maintenance and protection against end-to-end fusion of chromosomes. TRF2 is a component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection. TRF2 together with DCLRE1B/Apollo, is required to control the amount of DNA topoisomerase (TOP1, TOP2A and TOP2B), which is needed for telomere replication during fork passage and prevent aberrant telomere topology. TRF2 recruits TERF2IP/RAP1 to telomeres, thereby participating in to repressing homology-directed repair (HDR), which can affect telomere length.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage


***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87020-1-RR (TERF2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.