

For Research Use Only

MKRN2 Recombinant monoclonal antibody

Catalog Number:87026-1-RR



Basic Information

Catalog Number: 87026-1-RR	GenBank Accession Number: BC015715	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 23609	CloneNo.: 252088F2
Source: Rabbit	UNIPROT ID: Q9H000	Recommended Dilutions: WB: 1:2000-1:10000 IF/ICC: 1:200-1:800
Isotype: IgG	Full Name: makorin ring finger protein 2	
Immunogen Catalog Number: AG2881	Calculated MW: 416 aa, 47 kDa	
	Observed MW: 47 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: human	WB : HeLa cells, HepG2 cells, HEK-293T cells, K-562 cells IF/ICC : HepG2 cells,

Background Information

MKRN2(Probable E3 ubiquitin-protein ligase makorin-2) is also named as RNF62. It is an E3 ubiquitin ligase that catalyzes the covalent attachment of ubiquitin moieties onto substrate proteins. It may be playing a functional role in negatively regulating neurogenesis via PI3K/Akt signaling.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87026-1-RR (MKRN2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using MKRN2 antibody (87026-1-RR, Clone: 252088F2) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).