

For Research Use Only

TRPM2 Recombinant monoclonal antibody, PBS Only

Catalog Number:87036-1-PBS



Basic Information

Catalog Number:

87036-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC041570

GeneID (NCBI):

7226

UNIPROT ID:

O94759

Full Name:

transient receptor potential cation channel, subfamily M, member 2

Calculated MW:

171 kDa

Observed MW:

171 kDa

Purification Method:

Protein A purification

CloneNo.:

251122F2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MDA-MB-231 cells were subjected to SDS PAGE followed by western blot with 87036-1-RR (TRPM2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87036-1-PBS in a different storage buffer formulation.