

For Research Use Only

MYD88 Recombinant monoclonal antibody

Catalog Number: 87040-2-RR



Basic Information

Catalog Number: 87040-2-RR	GenBank Accession Number: BC013589	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 4615	CloneNo.: 252132B7
Source: Rabbit	UNIPROT ID: Q99836	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: myeloid differentiation primary response gene (88)	
Immunogen Catalog Number: AG19770	Calculated MW: 33 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human, rat	WB : Jurkat cells, MOLT-4 cells, K-562 cells, Raji cells, HepG2 cells, rat liver tissue

Background Information

MYD88 is a cytosolic adapter protein that plays a central role in the innate and adaptive immune response. MYD88 is the key adaptor protein for the interleukin-1 receptor (IL-1R) and most Toll-like receptors (TLRs), which are essential for innate immunity and pathogen-associated molecular pattern recognition. Mutations in the MyD88 gene lead to the development of cancer in humans and mice suggesting that MyD88 also plays a cell autonomous role in tissue homeostasis.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87040-2-RR (MYD88 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.