

For Research Use Only

MYD88 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87040-4-PBS



Basic Information

Catalog Number: 87040-4-PBS	GenBank Accession Number: BC013589	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4615	CloneNo.: 252132A6
Source: Rabbit	UNIPROT ID: Q99836	
Isotype: IgG	Full Name: myeloid differentiation primary response gene (88)	
Immunogen Catalog Number: AG19770	Calculated MW: 33 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Product Information

87040-4-PBS targets MYD88 as part of a matched antibody pair:

MP03184-1: 87040-4-PBS capture and 87040-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MYD88 concentration was determined to be 907.15 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load, 2,599.63 ng/mL in Raji cell extract based on a 1.20 mg/mL extract load and 123.62 ng/mL in HCT 116 cell extract based on a 1.30 mg/mL extract load.



Sandwich ELISA standard curve of MP03184-1, Human MYD88 Recombinant Matched Antibody Pair - PBS only. 87040-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19770. 87040-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL