

For Research Use Only

UBL3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87070-3-PBS



Basic Information

Catalog Number:

87070-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5932

GenBank Accession Number:

NM_007106.3

GeneID (NCBI):

5412

UNIPROT ID:

O95164

Full Name:

ubiquitin-like 3

Calculated MW:

13 kDa

Observed MW:

13 kDa

Purification Method:

Protein A purification

CloneNo.:

252163G4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Ubiquitin-like 3 (UBL3)/membrane-anchored Ub-fold protein (MUB) acts as a post-translational modification (PTM) factor to regulate efficient protein sorting to sEVs (small extracellular vesicles). (PMID: 31363817, PMID: 30258067).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87070-3-RR (UBL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87070-3-PBS in a different storage buffer formulation.