

For Research Use Only

CABYR Recombinant monoclonal antibody, PBS Only

Catalog Number:87074-1-PBS



Basic Information

Catalog Number:

87074-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3012

GenBank Accession Number:

BC011996

GeneID (NCBI):

26256

UNIPROT ID:

O75952

Full Name:

calcium binding tyrosine-(Y)-phosphorylation regulated

Calculated MW:

53 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

252259B4

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87074-1-RR (CABYR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87074-1-PBS in a different storage buffer formulation.

