

For Research Use Only

Rat CD45 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87083-1-PBS



Basic Information

Catalog Number: 87083-1-PBS	GenBank Accession Number: NM_138507.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 24699	CloneNo.: 252198C12
Source: Rabbit	UNIPROT ID: P04157-1	
Isotype: IgG	Full Name: protein tyrosine phosphatase, receptor type, C	
Immunogen Catalog Number: EG4666	Calculated MW: 143 kDa	
	Observed MW: 180 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
rat

Background Information

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87083-1-RR (CD45/PTPRC antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87083-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87083-1-RR against Rat CD45 were performed. The affinity constant is 14.7 nM.