

For Research Use Only

OSBPL10 Recombinant monoclonal antibody

Catalog Number: 87087-1-RR



Basic Information

Catalog Number: 87087-1-RR	GenBank Accession Number: BC003168	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 114884	CloneNo.: 252248A4
Source: Rabbit	UNIPROT ID: Q9BXB5	Recommended Dilutions: WB: 1:2000-1:10000 IF/ICC: 1:50-1:500
Isotype: IgG	Full Name: oxysterol binding protein-like 10	
Immunogen Catalog Number: AG7861	Calculated MW: 84 kDa	
	Observed MW: 75-80 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HuH-7 cells, PANC-1 cells, HepG2 cells, A431 cells, HT-29 cells, rat brain tissue
Species Specificity: human, rat	IF/ICC : HepG2 cells,

Background Information

OSBPL10, also known as ORP10, is a member of the oxysterol-binding protein family. OSBPL10 localizes on microtubules and acts as a regulator of cellular lipid biosyntheses and apolipoprotein B-100. Genetic variation in OSBPL10 is associated with dyslipidias and peripheral artery disease.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87087-1-RR (OSBPL10 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using OSBPL10 antibody (87087-1-RR, Clone: 252248A4) at dilution of 1:200 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002).