

For Research Use Only

TRPV3 Recombinant monoclonal antibody, PBS Only

Catalog Number:87098-1-PBS



Basic Information

Catalog Number:

87098-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC104868

GeneID (NCBI):

162514

UNIPROT ID:

Q8NET8

Full Name:

transient receptor potential cation channel, subfamily V, member 3

Calculated MW:

791 aa, 91 kDa

Observed MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

250990D11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Transient receptor potential cation channel, subfamily V, member 3, also known as TRPV3, is a human gene encoding the protein of the same name. The TRPV3 protein belongs to a family of nonselective cation channels that function in a variety of processes, including temperature sensation and vasoregulation. The thermosensitive members of this family are expressed in subsets of human sensory neurons that terminate in the skin, and are activated at distinct physiological temperatures. The TRPV3 channel has wide tissue expression that is especially high in the skin (keratinocytes) but also in the brain. It functions as a molecular sensor for innocuous warm temperatures.(PMID 12016205)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87098-1-RR (TRPV3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87098-1-PBS in a different storage buffer formulation.