

For Research Use Only

Phospho-TAU (Thr212/Ser214) Recombinant monoclonal antibody, PBS Only



Catalog Number: 87103-1-PBS

Basic Information

Catalog Number: 87103-1-PBS	GenBank Accession Number: BC000558	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4137	CloneNo.: 252171D2
Source: Rabbit	UNIPROT ID: P10636	
Isotype: IgG	Full Name: microtubule-associated protein tau	
	Calculated MW: 37-46, 79-81 kDa	
	Observed MW: 45-60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
human, mouse

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated mouse brain tissue, phosphatase inhibitor treated mouse brain tissue and λ phosphatase treated mouse brain tissue were subjected to SDS PAGE followed by western blot with 87103-1-RR (Phospho-TAU (Thr212/Ser214) antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Alpha Tubulin (66031-1-Ig) antibody as a loading control. This data was developed using the same antibody clone

