

For Research Use Only

USP3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87122-1-PBS



Basic Information

Catalog Number: 87122-1-PBS	GenBank Accession Number: BC018113	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9960	CloneNo.: 252249C3
Source: Rabbit	UNIPROT ID: Q9Y6I4	
Isotype: IgG	Full Name: ubiquitin specific peptidase 3	
Immunogen Catalog Number: AG14594	Calculated MW: 520 aa, 59 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

USP3(Ubiquitin carboxyl-terminal hydrolase 3) belongs to the peptidase C19 family. It may regulate the DNA damage response (DDR) checkpoint through deubiquitination of H2A at DNA damage sites. USP3 is implicated as a novel chromatin modifier in the maintenance of genome integrity. This protein has 2 isoforms produced by alternative splicing.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87122-1-RR (USP3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87122-1-PBS in a different storage buffer formulation.