

For Research Use Only

PDZD8 Recombinant monoclonal antibody

Catalog Number: 87134-1-RR



Basic Information

Catalog Number: 87134-1-RR	GenBank Accession Number: BC028375	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 118987	CloneNo.: 252321A11
Source: Rabbit	UNIPROT ID: Q8NEN9	Recommended Dilutions: WB: 1:1000-1:4000
Isotype: IgG	Full Name: PDZ domain containing 8	
Immunogen Catalog Number: AG21914	Calculated MW: 1154 aa, 129 kDa	
	Observed MW: 150-160 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HCT 116 cells, HeLa cells, HEK-293 cells
Species Specificity: human	

Background Information

PDZD8 (PDZ domain containing 8) is an SMP (synaptotagmin-like mitochondrial-lipid-binding) domain-containing protein belonging to the tubular lipid-binding protein (TULIP) superfamily of lipid transfer proteins (PMID: 37961523). PDZD8 is an integral ER transmembrane protein present at ER-mitochondria contacts (PMID: 29097544). It is involved in regulating ER-mitochondria associations, mitochondrial calcium ion homeostasis, and lipid transport between the ER and late endosomes/lysosomes (PMID: 37247681; 29097544; 33912962).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87134-1-RR (PDZD8 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87134-1-RR against Human PDZD8 were performed. The affinity constant is 0.398 nM.