

For Research Use Only

SAA1+SAA2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87185-1-PBS



Basic Information

Catalog Number:

87185-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4687

GenBank Accession Number:

GeneID (NCBI):

6288

UNIPROT ID:

PODJ18

Full Name:

serum amyloid A1

Calculated MW:

14kDa

Purification Method:

Protein A purification

CloneNo.:

251324C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

87185-1-PBS targets SAA1+SAA2 as part of a matched antibody pair:

MP02910-1: 87185-1-PBS capture and 87184-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Serum Amyloid A concentration was determined to be 5.6 ug/mL in human plasma.



Cytometric bead array standard curve of MP02910-1, Serum Amyloid A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87185-1-PBS. Detection antibody: 87184-1-PBS. Standard: Eg4687. Range: 0.625-80 ng/mL.