

For Research Use Only

SLC12A1/NKCC2 Recombinant monoclonal antibody

Catalog Number: 87210-1-RR



Basic Information

Catalog Number:

87210-1-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5325

GenBank Accession Number:

BC040138

GeneID (NCBI):

6557

UNIPROT ID:

Q13621

Full Name:

solute carrier family 12 (sodium/potassium/chloride transporters), member 1

Calculated MW:

121 kDa

Observed MW:

140-180 kDa

Purification Method:

Protein A purification

CloneNo.:

252363D10

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB: 37°C incubated mouse kidney tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated and 37°C incubated mouse kidney tissue were subjected to SDS PAGE followed by western blot with 87210-1-RR (SLC12A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.