

For Research Use Only

# MYO18A Recombinant monoclonal antibody

Catalog Number: 87218-4-RR



## Basic Information

Catalog Number:

87218-4-RR

Size:

100ul , Concentration: 1000 µg/ml by 399687

Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC039612

GeneID (NCBI):

399687

UNIPROT ID:

Q92614

Full Name:

myosin XVIII A

Calculated MW:

233 kDa

Observed MW:

230-245 kDa

Purification Method:

Protein A purification

CloneNo.:

251878A9

Recommended Dilutions:

WB: 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HeLa cells, HT-29 cells, K-562 cells, PC-12 cells

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87218-4-RR (MYO18A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.