

DNMT2 Recombinant monoclonal antibody

Catalog Number: 87236-1-RR

Basic Information

Catalog Number: 87236-1-RR	GenBank Accession Number: BC047733	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 1787	CloneNo.: 252459G1
Source: Rabbit	UNIPROT ID: O14717	Recommended Dilutions: WB: 1:1000-1:6000
Isotype: IgG	Full Name: tRNA aspartic acid methyltransferase 1	
Immunogen Catalog Number: AG5571	Calculated MW: 45 kDa	
	Observed MW: 43-45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human	WB : HCT 116 cells, HeLa cells, HEK-293 cells, K-562 cells

Background Information

Dnmt2 was initially assigned a member of the DNA methyltransferase family on the basis of its extensive homology with eukaryotic and prokaryotic DNA-(cytosine C5)-methyltransferases. However, in the apparent absence of a phenotype in dnmt2 knockout cells, Dnmt2's possible biological function remained unknown, even though Dnmt2 is strongly conserved and it is found in species ranging from *Schizosaccharomyces pombe* to human. Later very weak, residual DNA methylation activity was found with enzymes from different species. The finding that Dnmt2 is an active RNA methyltransferase capable of methylating the C38 position of the tRNA Asp came as a surprise. However, still no cellular function of the tRNA Asp methylation has been found, although in Zebrafish Dnmt2 knock-down caused a developmental phenotype. It is very intriguing that an enzyme that looks like a DNA methyltransferase can methylate RNA, in particular since the RNA and DNA specific m5C methyltransferases use different catalytic residues and a different mechanism for the methyl transfer reaction. The predicted molecular weight of Dnmt2 is 45 kDa, but it may migrate to the 40 kDa (PMID: 18567810)

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87236-1-RR (DNMT2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.