

For Research Use Only

Mouse CD316/IGSF8 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87310-2-PBS



Basic Information

Catalog Number: 87310-2-PBS	GenBank Accession Number: AF439263.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 140559	CloneNo.: 252500E12
Source: Rabbit	UNIPROT ID: Q8R366	
Isotype: IgG	Full Name: immunoglobulin superfamily, member 8	
Immunogen Catalog Number: EG6320	Calculated MW: 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

87310-2-PBS targets CD316/IGSF8 as part of a matched antibody pair:

MP02958-1: 87310-1-PBS capture and 87310-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IGSF8/CD316 concentration was determined to be 865.1 ng/mL in mouse brain tissue extract based on a 2.2 mg/mL extract load, A431 cell extract were served as a negative control.



Cytometric bead array standard curve of MP02958-1, MOUSE IGSF8/CD316 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87310-1-PBS. Detection antibody: 87310-2-PBS. Standard: Eg6320. Range: 0.625-80 ng/mL.