

For Research Use Only

Phospho-P53 (Ser33) Recombinant monoclonal antibody

Catalog Number:87312-1-RR



Basic Information

Catalog Number: 87312-1-RR	GenBank Accession Number: NML_000546	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7157	CloneNo.: 252563C6
Source: Rabbit	UNIPROT ID: P04637	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: tumor protein p53	
	Calculated MW: 44 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : etoposide treated T-47D cells,
Species Specificity: human	

Background Information

TP53, also known as P53 and NY-CO-13, belongs to the p53 family and has 9 isoforms. In SDS-Page, the observed molecular weight is about 53 kDa. Phosphorylation of p53 at specific amino acid residues plays a crucial role in mediating its various functions, including its role as a transcription factor that regulates the expression of genes involved in the cellular response to stress and DNA damage. Upon severe DNA damage, p53 is phosphorylated at specific amino acid residues, becomes stabilized, and activated to induce apoptosis-related genes. To date, 10 kinases have been identified to phosphorylate p53 for apoptosis at specific serine/threonine residues in the N- and C-terminus domains. (PMID: 30191640)

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Non-treated T-47D cells and etoposide treated T-47D cells were subjected to SDS PAGE followed by western blot with 87312-1-RR (Phospho-P53 (Ser33) antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with Lamin B1 (66095-1-Ig) antibody as a loading control.