

For Research Use Only

B3GALT6 Recombinant monoclonal antibody, PBS Only

Catalog Number:87370-1-PBS



Basic Information

Catalog Number:

87370-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_080605

GeneID (NCBI):

126792

UNIPROT ID:

Q96L58

Full Name:

UDP-Gal:betaGal beta 1,3-galactosyltransferase polypeptide 6

Calculated MW:

37 kDa

Observed MW:

37-40 kDa

Purification Method:

Protein A purification

CloneNo.:

252785F4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87370-1-RR (B3GALT6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87370-1-PBS in a different storage buffer formulation.