

For Research Use Only

B3GALT6 Recombinant monoclonal antibody

Catalog Number: 87370-1-RR



Basic Information

Catalog Number:

87370-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_080605

GeneID (NCBI):

126792

UNIPROT ID:

Q96L58

Full Name:

UDP-Gal:betaGal beta 1,3-galactosyltransferase polypeptide 6

Calculated MW:

37 kDa

Observed MW:

37-40 kDa

Purification Method:

Protein A purification

CloneNo.:

252785F4

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue, mouse liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87370-1-RR (B3GALT6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.