

For Research Use Only

OTX2 Recombinant monoclonal antibody

Catalog Number: 87372-1-RR



Basic Information

Catalog Number:

87372-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4323

GenBank Accession Number:

BC032579

GeneID (NCBI):

5015

UNIPROT ID:

P32243

Full Name:

orthodenticle homeobox 2

Calculated MW:

297 aa, 32 kDa

Observed MW:

31-35 kDa

Purification Method:

Protein A purification

CloneNo.:

252762H3

Recommended Dilutions:

WB: 1:5000-1:50000

IF/ICC: 1:250-1:1000

Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : ARPE-19 cells, mouse retina tissue, rat retina tissue

IF/ICC : Y79 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87372-1-RR (OTX2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed Y79 cells using OTX2 antibody (87372-1-RR, Clone: 252762H3) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Biolayer interferometry (BLI) kinetic assays of 87372-1-RR against Human OTX2 were performed. The affinity constant is 0.132 nM.