

For Research Use Only

CCL4L1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87376-1-PBS



Basic Information

Catalog Number: 87376-1-PBS	GenBank Accession Number: N/A	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9560	CloneNo.: 251685H1
Source: Rabbit	UNIPROT ID: Q8NHW4-1	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 4-like 1	
Immunogen Catalog Number: EG5335	Calculated MW: 10 kDa	
	Observed MW: 15 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PMA, LPS and Brefeldin A treated THP-1 cells were subjected to SDS PAGE followed by western blot with 87376-1-RR (CCL4L1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87376-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87376-1-RR against Human CCL4L1 were performed. The affinity constant is 4.87 nM.