

For Research Use Only

Human IgM Recombinant monoclonal antibody, PBS Only

Catalog Number: 87380-10-PBS



Basic Information

Catalog Number:

87380-10-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6742

GenBank Accession Number:

BC020240

GeneID (NCBI):

3507

Full Name:

immunoglobulin heavy constant mu

Purification Method:

Protein A purification

CloneNo.:

253289A9

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Specificity of the recombinant rabbit anti-human IgM antibody was confirmed by ELISA. The antibody specifically recognizes both human IgM monomer and pentamer, with no cross-reactivity observed against other human immunoglobulins (IgG, IgD, IgE) or IgM from mouse, pig, rat, and cat.