

For Research Use Only

Human IgM Recombinant monoclonal antibody

Catalog Number: 87380-10-RR



Basic Information

Catalog Number:

87380-10-RR

GenBank Accession Number:

BC020240

Purification Method:

Protein A purification

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

GeneID (NCBI):

3507

CloneNo.:

253289A9

Source:

Rabbit

Full Name:

immunoglobulin heavy constant mu

Isotype:

IgG

Immunogen Catalog Number:

EG6742

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Specificity of the recombinant rabbit anti-human IgM antibody was confirmed by ELISA. The antibody specifically recognizes both human IgM monomer and pentamer, with no cross-reactivity observed against other human immunoglobulins (IgG, IgD, IgE) or IgM from mouse, pig, rat, and cat.