

For Research Use Only

Human IgM Recombinant monoclonal antibody

Catalog Number: 87380-8-RR



Basic Information

Catalog Number: 87380-8-RR	GenBank Accession Number: BC020240	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3507	CloneNo.: 253289F4
Source: Rabbit	Full Name: immunoglobulin heavy constant mu	
Isotype: IgG		
Immunogen Catalog Number: EG6742		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Specificity of the recombinant rabbit anti-human IgM antibody was confirmed by ELISA, demonstrating exclusive reactivity to human IgM monomer with no cross-reactivity to other human immunoglobulins or IgM from mouse, pig, rat, or cat.