

For Research Use Only

# CD297/ART4 Recombinant monoclonal antibody

Catalog Number: 87399-1-RR



## Basic Information

<b>Catalog Number:</b> 87399-1-RR	<b>GenBank Accession Number:</b> NM_026639.2	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 109978	<b>CloneNo.:</b> 252721F4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9CRA0	<b>Recommended Dilutions:</b> WB: 1:1000-1:4000
<b>Isotype:</b> IgG	<b>Full Name:</b> ADP-ribosyltransferase 4	
<b>Immunogen Catalog Number:</b> EG6315	<b>Calculated MW:</b> 33 kDa	
	<b>Observed MW:</b> 33 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse spleen tissue, mouse thymus tissue, mouse liver tissue
<b>Species Specificity:</b> mouse	

## Background Information

CD297/ART4 is a GPI-anchored ecto-ADP-ribosyltransferase best known as the molecular carrier of Dombrock blood group antigens on erythrocytes. Expressed on red blood cells, hematopoietic progenitors, and select immune cells, it participates in cell-surface signaling and immune recognition, though its catalytic activity is limited. CD297/ART4 polymorphisms underlie clinically significant transfusion reactions, making it a key target in hematology, transfusion medicine, and erythroid biology.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

