

CD41 Recombinant monoclonal antibody

Catalog Number: 87447-4-RR

Basic Information

Catalog Number: 87447-4-RR	GenBank Accession Number: 3674	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 3674	CloneNo.: 252680C11
Source: Rabbit	UNIPROT ID: P08514	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: integrin, alpha 2b (platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41)	
Immunogen Catalog Number: EG10496	Observed MW: 110-125 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human peripheral blood platelets,
Species Specificity: human	

Background Information

The integrins are a superfamily of cell adhesion receptors that bind to extracellular matrix ligands, cell-surface ligands, and soluble ligands (PMID: 17543136). They are transmembrane $\alpha\beta$ heterodimers and at least 24 distinct integrin heterodimers are formed by the combination of 18 α and eight β known subunits (PMID: 17543136; 20029421). In addition to mediating cell adhesion, integrins also play important roles in modulating signal transduction pathways that control cellular responses, including migration, proliferation, differentiation, and apoptosis (PMID: 19118207). Integrin α -IIb (ITGA2B, CD41) is expressed on platelets, megakaryocytes and some hematopoietic progenitor cells (PMID: 11934866). This protein undergoes post-translational cleavage to yield disulfide linked light beta (25 kDa) and heavy alpha (125 kDa) chains, which join with integrin beta3 (CD61) to form a receptor for fibronectin, fibrinogen, plasminogen, prothrombin, thrombospondin and vitronectin. CD41/CD61 is crucial for platelet aggregation through binding of soluble fibrinogen. Jurkat and U937 cell lines are negative controls for CD41 (PMID: 2351656).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87447-4-RR (ITA2B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.