

For Research Use Only

TXNIP Recombinant monoclonal antibody, PBS Only

Catalog Number: 87483-1-PBS



Basic Information

Catalog Number: 87483-1-PBS	GenBank Accession Number: BC093702	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10628	CloneNo.: 252846D1
Source: Rabbit	UNIPROT ID: Q9H3M7	
Isotype: IgG	Full Name: thioredoxin interacting protein	
Immunogen Catalog Number: AG33838	Calculated MW: 391 aa, 44 kDa	
	Observed MW: 44-55 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

TXNIP, also known as VDUP-1 or TBP-2, belongs to alpha-arrestin protein family and is perhaps the only family member known to bind thioredoxin (TRX). TXNIP was induced by Vitamin D3, but not induced by another monocyte or macrophage differentiation inducer: phorbol 12-myristate 13-acetate (PMA). TXNIP bound catalytic active center of thioredoxin (TRX), which protected cells against oxidative stress. TXNIP was found to be a negative regulator of thioredoxin activity and inducer of the intracellular level of reactive oxygen species (ROS). TXNIP plays an important role in a wide variety of biological functions, such as the regulation of cell death, growth, differentiation, and energy metabolism.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87483-1-RR (TXNIP antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87483-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87483-1-RR against Human TXNIP were performed. The affinity constant is below 1 pM.