

For Research Use Only

TUBG2 Recombinant monoclonal antibody

Catalog Number: 87516-1-RR



Basic Information

Catalog Number: 87516-1-RR	GenBank Accession Number: BC051890	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 27175	CloneNo.: 252833C10
Source: Rabbit	UNIPROT ID: Q9NRH3	Recommended Dilutions: WB: 1:2000-1:10000 IF/ICC: 1:250-1:1000
Isotype: IgG	Full Name: tubulin, gamma 2	
Immunogen Catalog Number: AG27728	Calculated MW: 451 aa, 51 kDa	
	Observed MW: 51 kDa	

Applications

Tested Applications:
WB, IF/ICC, ELISA

Species Specificity:
human, mouse, rat, canine

Positive Controls:

WB : HeLa cells, T-47D cells, SK-BR-3 cells, MDCK cells, COS-7 cells, mouse cerebellum tissue, rat cerebellum tissue

IF/ICC : HeLa cells,

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87516-1-RR (TUBG2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using TUBG2 antibody (87516-1-RR, Clone: 252833C10) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).