

For Research Use Only

TBP Recombinant monoclonal antibody

Catalog Number: 87550-1-RR



Basic Information

Catalog Number: 87550-1-RR	GenBank Accession Number: BC110341	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6908	CloneNo.: 252942E3
Source: Rabbit	UNIPROT ID: P20226	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:500-1:2000
Isotype: IgG	Full Name: TATA box binding protein	
Immunogen Catalog Number: AG12383	Calculated MW: 338 aa, 38 kDa	
	Observed MW: 37-45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human	WB : HEK-293 cells, HeLa cells, PC-3 cells, A2780 cells, COLO 320 cells
	IF/ICC : HepG2 cells,

Background Information

TBP, TATA box-binding protein, is a DNA-binding subunit of TFIID, which is a multisubunit complex essential for the expression of most protein-encoding genes. The well-conserved C-terminal domains of 180 aa is both necessary and sufficient for interaction with DNA and for assembly of the basal transcription apparatus. The N terminus of TBP modulates the DNA-binding activity of the C terminus of the protein. Once binds to the TATA box mediated by TBP, TFIID initiate the transcription step of the pre-initiation complex (PIC), playing a role in the activation of eukaryotic genes transcribed by RNA polymerase II.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87550-1-RR (TBP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using TBP antibody (87550-1-RR, Clone: 252942E3) at dilution of 1:1000 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).