

For Research Use Only

SLC25A29 Recombinant monoclonal antibody

Catalog Number: 87551-1-RR



Basic Information

Catalog Number: 87551-1-RR	GenBank Accession Number: BC132966	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 123096	CloneNo.: 252932H6
Source: Rabbit	UNIPROT ID: Q8N8R3	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Full Name: solute carrier family 25, member 29	
Immunogen Catalog Number: AG24662	Calculated MW: 303 aa, 32 kDa	
	Observed MW: 32 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87551-1-RR (SLC25A29 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.