

For Research Use Only

CNTN4 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87572-1-PBS



Basic Information

Catalog Number: 87572-1-PBS	GenBank Accession Number: NM_001206955.2	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 152330	CloneNo.: 252031B8
Source: Rabbit	UNIPROT ID: Q8IWW2-1	
Isotype: IgG	Full Name: contactin 4	
Immunogen Catalog Number: EG5574	Calculated MW: 113 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Contactins are nervous system-specific members of the immunoglobulin superfamily, comprising six subtypes: contactin, TAG-1, BIG-1, BIG-2, NB-2 and NB-3 (PMID: 19262165). They are structurally and functionally conserved among species and predominantly expressed in neurons (PMID: 25300137). As a GPI-anchored synaptic adhesion molecule, CNTN4 (BIG-2) regulates brain axonal activities (PMID: 14571131). Elevated CNTN4 expression in tumors suppresses the anti-tumor immunity of T cells (PMID: 39392894).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 87572-1-RR (CNTN4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87572-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87572-1-RR against Human CNTN4 were performed. The affinity constant is below 1 pM.