

For Research Use Only

TMEM150 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87578-1-PBS



Basic Information

Catalog Number:

87578-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36287

GenBank Accession Number:

BC050466

GeneID (NCBI):

129303

UNIPROT ID:

Q86TG1

Full Name:

transmembrane protein 150

Calculated MW:

271 aa, 29 kDa

Observed MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

252959H6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87578-1-RR (TMEM150 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87578-1-PBS in a different storage buffer formulation.